

The kit contains selected basic components and special aids that are designed to be distributed to people seriously affected by emergency situations (floods, earthquakes, war conflicts, etc.).

The kit provides immediate emergency relief, without need of further assistance.

In addition, the humanitarian aid kit includes 0.5 l of drinking water ready to use.

Individual components of the kit are stored in a portable plastic container with a safety lid, which is designed as a measuring container that can be also used for preparation of drinking water.



MEDICAL KIT: Panadol – pain relief 10 tbl., Loperon – diarrhoea relief 10 tbl., Ajatin 50 ml – disinfectant 1 pc, Adhesive plaster – divided 6 x 2 cm 10 pcs, Sterile gauze 7.5 x 7.5 cm 2 pcs, Triangular bandage 1 pc, Non-sterile gloves 1 pair, Safety pin 2 pcs, Klickbox – plastic box 1 pc

WATER TREATMENT KIT: Aqua steril – drinking water treatment chemical - 1 set, Plastic box with a lid and a handle (volume 5.6 l) - 1 pc

Water that has been treated using this kit is suitable for drinking and cooking; it may be used to prepare drinks, wash food, clean wounds, as well as for personal hygiene, etc.

The kit enables you to prepare safe drinking water and especially to remove bacteria, viruses and parasites contained in the raw water in places where the quality of water sources is unknown or hygienically unsafe.

The treated water remains conserved for a long time. The kit disinfects raw water of unknown or poor quality and can be used instead of boiling the water in order to ensure microbial safety. Using treated water prevents gastrointestinal conditions that may be caused due to use of unsafe water

PERSONAL SAFETY KIT: Isothermic emergency blanket 1 pc, Drinking water 125 ml 4 pcs, Raincoat 1 pc, Pencil 1 pc, Note pad 1 pc, Whistle 1 pc, String 1 pc, Non-sterile gloves 1 pair, Cup candle 3 pcs

COOKING SET: Solid fuel cooker 1 pc, Solid fuel – tablets 12 pcs, Matches 1 box, Cup 1 pc, Plastic cutlery 2 pcs, Non-sterile gloves 1 pair, PE bag 1 pc, PE clip bag 5 pcs, Bag for contaminated food 2 pcs, Storage: in a dry place, at temperatures 15 – 25 °C

HYGIENIC SET: Antibacterial soap 1pc, Paper tissue 4 packets, Antibacterial tissue 3 pcs, Rubber band 3 pcs

Effectiveness of the new hygienic antibacterial tissues: Nanosilver technology enables effective removal of bacteria and yeasts. Effective for 12 hours! Materials treated with the nanosilver technology acquire new antibacterial properties and ensure reliable and permanent removal of bacteria, yeasts and moulds, e.g. Staphylococcus aureus (MRSA), beta haemolytic streptococcus, enterococcus, Listeria sp., Salmonella sp., Escherichia coli, Klebsiella sp., Enterobacter cloacae, Proteus sp., Morganella morganii, Providencia stuartii, Citrobacter sp., Pseudomonas aeruginosa, Stenotrophomonas gonorrhoeae, Chlamydia trachomatis, Bacillus subtilis, Candida albicans.